

## Part I

# Materials and Methods

## 1 Buffers

### 1.1 Growth Media

Prepare required volume of Luria Broth. Add appropriate antibiotics to following concentration

| Antibiotic      | Concentration ( $\mu\text{g ml}^{-1}$ ) |
|-----------------|---|
| Ampicillin      | 100                                     |
| Chloramphenicol | 39                                      |

### 1.2 Induction Media

Prepare growth media with appropriate antibiotic. Add appropriate inducers to following concentration

| Inducer   | Concentration              |
|-----------|----------------------------|
| Arabinose | $100 \mu\text{g ml}^{-1}$  |
| IPTG      | $25 \mu\text{mol dm}^{-3}$ |

### 1.3 5X Isothermal Buffer

|                     | Stock Concentration (mmol) | Volume ( $\mu$ l) |
|---------------------|----------------------------|-------------------|
| PEG-8000            | 25%                        | 0.75 g            |
| Tris-HCl pH 7.5     | 500                        | 1500              |
| MgCl <sub>2</sub>   | 50                         | 75                |
| DTT                 | 50                         | 150               |
| dATP                | 1                          | 30                |
| dTTP                | 1                          | 30                |
| dCTP                | 1                          | 30                |
| dGTP                | 1                          | 30                |
| NAD                 | 5                          | 300               |
| Nuclease free water |                            | ...               |
| <b>Total</b>        | 5X                         | 3000              |

Use nuclease free water as necessary to make up to 3 ml.

### 1.4 Gibson Master Mix

Prepare on ice in the ligase tube.

|                      | Stock Concentration (U $\mu$ l <sup>-1</sup> ) | Volume ( $\mu$ l) |
|----------------------|--|-------------------|
| Taq ligase           | 40   | 50                |
| 5X Isothermal Buffer | 5X   | 100               |
| T5 Exonuclease       | 1  | 2                 |
| Phusion Polymerase   | 2  | 6.25              |
| Nuclease free water  |  | 216.75            |
| <b>Total</b>         | 1.33X  | 375               |

For best performance, aliquot 25  $\mu$ l to 75  $\mu$ l portions in 1.5 ml microcentrifuge tubes and store at  $-20^{\circ}$ .

## 2 Cell Protocols

### 2.1 Overnight Cultures

Prepare in a 14ml shaky top tube 2 ml growth media with appropriate antibiotic. Inoculate a single colony in the growth media and grow overnight at 37°.

### 2.2 Inducing Cultures

Prepare 2 µl overnight culture. Prepare in a 14 µl shaky top tube 2 ml induction media with appropriate inducer. Inoculate induction media with 20 µl overnight culture.

### 2.3 Preparation of CaCl Competent Cells

Prepare 2 ml overnight culture. Meanwhile, prewarm 500 ml LB and two 1 l or one 2 l conical flasks and pre-cool a centrifuge to 4°. Also cool 16 1.5 ml microcentrifuge tubes. Add the overnight culture to the warmed 500 ml in the conical flask(s) and incubate in a shaking incubator for 2 h to 3 h. Incubate cells on ice for 10 min. Spin down cells for 3 min at 6krpm. Resuspend cells in 10 ml 0.1 mmol dm<sup>-3</sup> CaCl. Incubate cells on ice for 10 min. Spin down cells for 3 min at 6krpm. Resuspend cells in 5 ml 0.1 mmol dm<sup>-3</sup> CaCl in 15% glycerol. Aliquot 300 µl into each of the cooled microcentrifuge tubes. Freeze at -80°.

### 2.4 Transformation of CaCl Competent Cells

Thaw 20 µl to 200 µl competent cells on ice in 1.5 ml microcentrifuge tubes. Add 1 µl plasmid DNA or 20 µl ligation assembly DNA or complete Gibson Assembly DNA. Incubate on ice for 30 min. Heat shock for 45 s to 120 s at 42° Incubate on ice for at least 2 min. Add 800 µl LB. Incubate for 1 h at 37° in a shaking incubator. If required, spin down cells and remove excess LB to concentrate. Plate out 10 µl to 100 µl on LB agar plates with appropriate antibiotic. Incubate plates overnight at 37°.

## 2.5 Miniprep

## 2.6 Glycerol Stock

Prepare 2 ml overnight culture with appropriate antibiotic. Add 250  $\mu$ l to a 1.5 ml cryogenic tube. Add 1.15  $\mu$ l overnight culture and mix well. Freeze in  $-80^{\circ}$  freezer.

Source <http://web.wi.mit.edu/sive/pub/Lab>

# 3 DNA Protocols

## 3.1 Restriction Digest

Prepare on ice in order the following mixture.

|               | Volume( $\mu$ l) |
|---------------|------------------|
| DNA           | 50               |
| Buffer        | 6                |
| Milli Q Water | 1.4              |
| Enzyme A      | 1                |
| Enzyme B      | 1                |
| <b>Total</b>  | <b>60</b>        |

Incubate the mixture at  $37^{\circ}$  in a shaking incubator for 2-3 hours. If necessary, dephosphorylate the digested product.

## 3.2 Dephosphorylation

Prepare on ice in order the following mixture.

|               | Volume( $\mu$ l) |
|---------------|------------------|
| DNA           | 15               |
| SAP Buffer    | 2                |
| Milli Q Water | 2                |
| SAP           | 1                |
| <b>Total</b>  | <b>60</b>        |

Incubate the mixture at 37° for 1-2 hours.

### 3.3 PCR

### 3.4 Gibson Assembly

### 3.5 DNA Precipitation